

**Introducing..**

## **FloGripp™**

### **Flotec's New Portable Oxygen Cylinder Post-Valve Hand-Grip!**

Flotec's new **FloGripp™** replaces the sometimes burdensome "Toggle" lever on existing post valves which can easily allow a cylinder to inadvertently open and possibly discharge when stored, or during transport of the oxygen cylinders within soft-sided O<sub>2</sub> jump bags.



**FloGripp™** also eliminates the metallic securing pin used to hold the toggle in place That can cause cuts and abrasions during everyday use.

- Black anodized non-glare finish
- Knurled finish for enhanced **Gripp**-ability
- Attachment via twin side set screws with supplied allen wrench
- All Metal Construction using 6061 Heat-Tempered Aluminum Alloy
- Same outer diameter as post valve body for ease of regulator attachment
- Portable oxygen cylinder use has never been easier or more convenient!





## Post Valve Selection Chart

Complete Oxygen Cylinder Post Valve Assemblies with **FloGripp**

### For Aluminum Cylinders

PV830B	each	CGA 870 Valve (3/4-16 UNF), Knurled Aluminum Knob
PV831B	each	CGA 870 Valve (3/4-16 UNF), Knurled Aluminum Knob, LH Gauge
PV832B	each	CGA 870 Valve (3/4-16 UNF), Knurled Aluminum Knob, RH Gauge

### For Steel Cylinders

PV880B	each	CGA 870 Valve (1/2-14 NGT), Knurled Aluminum Knob
PV881B	each	CGA 870 Valve (1/2-14 NGT), Knurled Aluminum Knob, LH Gauge
PV882B	each	CGA 870 Valve (1/2-14 NGT), Knurled Aluminum Knob, RH Gauge



## Replacement Knob Assembly

*For conversion of existing toggle lever style oxygen cylinder post valve stems*

PVX30B	each	FloGripp Knob Assembly Kit
PVX30B-1	each	FloGripp Teflon Seating Washer
PVX30B-2	each	FloGripp Replacement Allen Wrench
PVX30B-3	each	Set Screws ( long )
PVX30B-4	each	Set Screws ( short )
PVX30B-5	each	Knob ( anodized )

