

----- HOSE ASSEMBLY LENGTH (Over all length (PER INCH))

----- FITTING END # 2

Select on two digit code from the group below

----- FITTING END # 1

Select on two digit code from the group below

3Y - DISS F Knurled Nut & Nipple w/ O-Ring Carbon Dioxide (grey)  
 3T - DISS F Knurled Nut & Nipple w/ O-Ring Medical Air (yellow)  
 3Q - DISS F Knurled Nut & Nipple w/ O-Ring Nitrogen (black)  
 3G - DISS F Knurled Nut & Nipple w/ O-Ring Nitrous Oxide (blue)  
 3S - DISS F Knurled Nut & Nipple w/ O-Ring Oxygen (green)  
 3C - DISS F Knurled Nut & Nipple w/ O-Ring Oxygen (white)  
 3W - DISS F Knurled Nut & Nipple w/ O-Ring VAC/Suction (white)  
 3E - DISS F Knurled Nut & Nipple w/ O-Ring WAGD/EVAC (purple)

4Y - DISS F Hex Nut & Nipple Carbon Dioxide (grey)  
 4U - DISS F Hex Nut & Nipple He-o2  
 4Z - DISS F Hex Nut & Nipple Instrument Air  
 4T - DISS F Hex Nut & Nipple Medical Air (yellow)  
 4Q - DISS F Hex Nut & Nipple Nitrogen (black)  
 4G - DISS F Hex Nut & Nipple Nitrous Oxide (blue)  
 4X - DISS F Hex Nut & Nipple o2-He  
 4S - DISS F Hex Nut & Nipple Oxygen (green)  
 4W - DISS F Hex Nut & Nipple VAC/Suction (white)  
 4E - DISS F Hex Nut & Nipple WAGD/EVAC (purple)

5Y - DISS F Hex Nut & Nipple w/ O-ring Carbon Dioxide (grey)  
 5U - DISS F Hex Nut & Nipple w/O-ring He-O2  
 5Z - DISS F Hex Nut & Nipple w/O-ring Instrument Air  
 5T - DISS F Hex Nut & Nipple w/O-ring Medical Air (yellow)  
 5Q - DISS F Hex Nut & Nipple w/O-ring Nitrogen (black)  
 5G - DISS F Hex Nut & Nipple w/O-ring Nitrous Oxide (blue)  
 5X - DISS F Hex Nut & Nipple w/O-ring o2-He  
 5S - DISS F Hex Nut & Nipple w/O-ring Oxygen (green)  
 5W - DISS F Hex Nut & Nipple w/O-ring VAC/Suction (white)  
 5E - DISS F Hex Nut & Nipple w/O-ring WAGD/EVAC (purple)

6S - DISS F Wing Nut & Nipple Oxygen (green)

7Y - DISS Male Body Adaptor w/ hand grip collar Carbon Dioxide (grey)  
 7T - DISS Male Body Adaptor w/ hand grip collar Medical Air (yellow)  
 7Q - DISS Male Body Adaptor w/ hand grip collar Nitrogen (black)  
 7G - DISS Male Body Adaptor w/ hand grip collar Nitrous Oxide (blue)  
 7S - DISS Male Body Adaptor w/ hand grip collar Oxygen (green)  
 7W - DISS Male Body Adaptor w/ hand grip collar VAC/Suction (white)  
 7E - DISS Male Body Adaptor w/ hand grip collar WAGD/EVAC (purple)

8Y - DISS Male w/ Demand Check, Nipple opens valve Carbon Dioxide (grey)  
 8T - DISS Male w/ Demand Check, Nipple opens valve Medical Air (yellow)  
 8Q - DISS Male w/ Demand Check, Nipple opens valve Nitrogen (black)  
 8G - DISS Male w/ Demand Check, Nipple opens valve Nitrous Oxide (blue)  
 8S - DISS Male w/ Demand Check, Nipple opens valve Oxygen (green)  
 8W - DISS Male w/ Demand Check, Nipple opens valve VAC/Suction (white)  
 8E - DISS Male w/ Demand Check, Nipple opens valve WAGD/EVAC (purple)

9T - DISS Male Body Adaptor Air  
 9U - DISS Male Body Adaptor He-o2  
 9Z - DISS Male Body Adaptor Instrument Air  
 9X - DISS Male Body Adaptor o2-He  
 9S - DISS Male Body Adaptor Oxygen

T3 - Schrader Style Male Adaptor w/ Swivel Stem Air  
 Q3 - Schrader Style Male Adaptor w/ Swivel Stem Nitrogen  
 G3 - Schrader Style Male Adaptor w/ Swivel Stem Nitrous Oxide  
 S3 - Schrader Style Male Adaptor w/ Swivel Stem Oxygen  
 W3 - Schrader Style Male Adaptor w/ Swivel Stem VAC/Suction

T2 - Schrader Style Female Coupler Medical Air  
 Q2 - Schrader Style Female Coupler Nitrogen  
 G2 - Schrader Style Female Coupler Nitrous Oxide  
 S2 - Schrader Style Female Coupler Oxygen  
 W2 - Schrader Style Female Coupler VAC/Suction

T1 - Schrader Style Male Adaptor non-swivel stem Medical Air  
 Q1 - Schrader Style Male Adaptor non-swivel stem Nitrogen  
 G1 - Schrader Style Male Adaptor non-swivel stem Nitrous Oxide  
 S1 - Schrader Style Male Adaptor non-swivel stem Oxygen  
 W1 - Schrader Style Male Adaptor non-swivel stem VAC/Suction

D1 - Chemetron Style Male Adaptor Medical Air (yellow)  
 N1 - Chemetron Style Male Adaptor Nitrous Oxide (blue)  
 C1 - Chemetron Style Male Adaptor Oxygen (green)  
 V1 - Chemetron Style Male Adaptor VAC/Suction (white)  
 E1 - Chemetron Style Male Adaptor WAGD/EVAC (purple)

D2 - Chemetron Style Female Coupler Medical Air (yellow)  
 N2 - Chemetron Style Female Coupler Nitrous Oxide (blue)  
 C2 - Chemetron Style Female Coupler Oxygen (green)  
 V2 - Chemetron Style Female Coupler VAC/Suction (white)  
 E2 - Chemetron Style Female Coupler WAGD/EVAC (purple)

R1 - Puritan Style Male Adaptor with HT collar Medical Air (yellow)  
 K1 - Puritan Style Male Adaptor with HT collar Nitrous Oxide (blue)  
 P1 - Puritan Style Male Adaptor with HT collar Oxygen (green)  
 M1 - Puritan Style Male Adaptor with HT collar VAC/Suction (white)  
 B1 - Puritan Style Male Adaptor with HT collar WAGD/EVAC (purple)

R2 - Puritan Style Female Coupler Medical Air (yellow)  
 K2 - Puritan Style Female Coupler Nitrous Oxide (blue)  
 P2 - Puritan Style Female Coupler Oxygen (green)  
 M2 - Puritan Style Female Coupler VAC/Suction (white)  
 B2 - Puritan Style Female Coupler WAGD/EVAC (purple)

Y1 - Ohmeda Style Male Adaptor Carbon Dioxide (grey)  
 A1 - Ohmeda Style Male Adaptor Medical Air (yellow)  
 H1 - Ohmeda Style Male Adaptor Nitrous Oxide (blue)  
 O1 - Ohmeda Style Male Adaptor Oxygen (green)  
 J1 - Ohmeda Style Male Adaptor VAC/Suction (white)  
 F1 - Ohmeda Style Male Adaptor WAGD/EVAC (purple)

Y2 - Ohmeda Style Female Coupler Carbon Dioxide (grey)  
 A2 - Ohmeda Style Female Coupler Medical Air (yellow)  
 H2 - Ohmeda Style Female Coupler Nitrous Oxide (blue)  
 O2 - Ohmeda Style Female Coupler Oxygen (green)  
 J2 - Ohmeda Style Female Coupler VAC/Suction (white)  
 F2 - Ohmeda Style Female Coupler WAGD/EVAC (purple)

1H - Male Hansen QD  
 2H - Female Hansen QD

22 - FEMALE 1/8" NPT PIPE THREAD  
 24 - FEMALE 1/4" NPT PIPE THREAD

12 - MALE 1/8" NPT PIPE THREAD  
 14 - MALE 1/4" NPT PIPE THREAD

00-NO END CONNECTION BLANK END

----- HOSE COLOR

B1 - Black (Nitrogen) Non-Conductive 1/4" ID Hose  
 B2 - Black (Nitrogen) Conductive 1/4" ID Hose  
 C1 - Clear w/ White Braid (General Purpose) Non-Conductive 1/4" ID Hose  
 G1 - Green (Oxygen) Non-Conductive 1/4" ID Hose  
 G2 - Green (Oxygen) Conductive 1/4" ID Hose  
 G3 - Green (Oxygen) Rubber 1/4" ID Hose  
 P1 - Purple (WAGD (EVAC)) Non-Conductive 1/4" ID Hose  
 P2 - Purple (WAGD (EVAC)) Conductive 1/4" ID Hose  
 R2 - Grey (Carbon Dioxide) Conductive 1/4" ID Hose  
 U1 - Blue (Nitrous Oxide) Non-Conductive 1/4" ID Hose  
 U2 - Blue (Nitrous Oxide) Conductive 1/4" ID Hose  
 W1 - White (Vacuum (Suction)) Non-Conductive 1/4" ID Hose  
 W2 - White (Vacuum (Suction)) Conductive 1/4" ID Hose  
 Y1 - Yellow (Medical Air) Non-Conductive 1/4" ID Hose  
 Y2 - Yellow (Medical Air) Conductive 1/4" ID Hose

----- PRODUCT GROUP:

HA - Hose Assembly