

FLOTEC INC

7625 West New York Street
Indianapolis, Indiana 46214-4911

Web Site: <http://www.floteco2.com>

E-mail: floteco@floteco2.com

Phone: 317-273-6960

Fax: 317-273-6979

Order desk: 1-800-401-1723

Order desk fax: 1-800-515-9254

Oxygen Adapter

Flotec's new **Oxygen Adapter** eliminates the problems of incompatible connections associated with different types of oxygen delivery systems allowing for hassle-free interchangeability between the two most common oxygen tubing connections utilized in medical facilities everyday.

Accepting both 22mm anesthesia / ventilator tubing as well as standard oxygen cannulas, the new **Oxygen Adapter**



Saves time, reduces waste and provides for a seamless changeover from one type of tubing connection to the other.

Patient transports, such as from surgery to recovery, usually require using not one but up to three different oxygen delivery tube / cannula set-ups during this process. With the **Oxygen Adapter** only one is needed.



Seamless continuous therapy between anesthesia, transport oxygen e-cylinder regulator and recovery room ventilator and /or wall mounted oxygen flow meter.

The **Oxygen Adapter** can be ordered factory installed as part of your next Flotec regulator order in lieu of the standard hose barb or DISS-1240 outlet or as a separate accessory threading onto the devices' male DISS-1240 outlet.

Part #	Description
INI-UOA-100	Male 1/8" NPT Threaded Connection
INI-UOA-200	Female DISS-1240 Threaded Connection