

Medical Oxygen Analyzers, Monitors & Sensors-ANI



[ANI AII 2000A Oxygen Analyzer](#)

Auto diagnostics, menu driven state-of-the-art electronics, large backlit LCD for easy viewing, integral tripod stand and dove tail bracket. Proprietary advanced galvanic sensor capable of oxygen analysis from 0-100% and 60 months expected life in air. General purpose analyzer operates 1,500 hours on 2 AA Alkaline batteries. Optional carrying case. Ranges - Standard: 0-100% Application - Spot checking respiratory therapy ventilators, neo-natal incubators, anesthesia



[ANI AII 2000HC Oxygen Analyzer](#)

Auto diagnostics, menu driven state-of-the-art electronics, internal oxygen sensor, large backlit LCD for easy viewing, integral tripod stand and dove tail bracket. Proprietary advanced galvanic sensor capable of oxygen analysis from 0-100% and 60 months expected life in air. General purpose analyzer operates 1,500 hours on 2 AA Alkaline batteries. Optional carrying case. Ranges - Standard: 0-100% Application - Spot checking tanks, concentrators and oxygen systems



[ANI-AII 2000M Oxygen Monitor](#)

Auto diagnostics, menu driven state-of-the-art electronics, audible and visual HI-LO alarms, LO alarm 15-99% oxygen, HI alarm 16-100% oxygen, large backlit LCD for easy viewing, integral tripod stand and dove tail bracket. Proprietary advanced galvanic sensor capable of oxygen analysis from 0-100% and 60 months expected life in air. General purpose monitor operates 1,500 hours on 2 AA Alkaline batteries. Optional carrying case. Ranges - Standard: 0-100% Application - Continuous use or spot check respiratory therapy ventilators, neo-natal incubators, anesthesia

<p align="center"><i>ANI-33</i></p>	<p align="center"><i>ANI-33-1</i></p>	<p align="center"><i>ANI-33-2</i></p>	<p align="center"><i>ANI-55</i></p>
			
<p>Range: 0-100% Oxygen Signal Output: 165-235 μA Response: 30 sec Expected Life: 24 months Warranty: 12 months</p>	<p>Range: 0-100% Oxygen Signal Output: 165-235 μA Response: 30 sec Expected Life: 24 months Warranty: 12 months</p>	<p>Range: 0-100% Oxygen Signal Output: 165-235 μA Response: 30 sec Expected Life: 24 months Warranty: 12 months</p>	<p>Range: 0-100% Oxygen Signal Output: 25-40 mV Response: 9 sec Expected Life: 42 months Warranty: 14 months</p>
<p align="center"><i>ANI-58</i></p>	<p align="center"><i>ANI-58HC</i></p>	<p align="center"><i>ANI-60</i></p>	<p align="center"><i>ANI-60-08</i></p>
			
<p>Range: 0-100% Oxygen Signal Output: 70-160 mV Response: 13 sec Expected Life: 60 months Warranty: 18 months</p>	<p>Range: 0-100% Oxygen Signal Output: 70-160 mV Response: 20 sec Expected Life: 60 months Warranty: 18 months</p>	<p>Range: 0-100% Oxygen Signal Output: Response:13 sec Expected Life: 60 months Warranty: 18 months</p>	<p>Range: 0-100% Oxygen Signal Output: Response:6 sec Expected Life: 60 months Warranty: 18 months</p>
<p align="center"><i>ANI-60-23</i></p>	<p align="center"><i>ANI-60-3</i></p>	<p align="center"><i>ANI-60HC</i></p>	<p align="center"><i>ANI-75-KE1</i></p>
			
<p>Range: 0-100% Oxygen Signal Output: Response:6 sec Expected Life: 60 months Warranty: 18 months</p>	<p>Range: 0-100% Oxygen Signal Output: Response:13 sec Expected Life: 60 months Warranty: 14 months</p>	<p>Range: 0-100% Oxygen Signal Output: Response:13 sec Expected Life: 60 months Warranty: 18 months</p>	<p>Range: 0-100% Oxygen Signal Output: 10-14 mV Response 90%: 13 sec Expected Life: 60 months Warranty: 18 months</p>

<i>ANI-75-KE10+</i>	<i>ANI-75-KE2</i>	<i>ANI-75-KE4</i>	<i>ANI-75-KE6</i>
			
Range: 0-100% Oxygen Signal Output: 10-14 mV Response 90%: 6 sec Expected Life: 60 months Warranty: 18 months	Range: 0-100% Oxygen Signal Output: 10-14 mV Response 90%: 13 sec Expected Life: 60 months Warranty: 18 months	Range: 0-100% Oxygen Signal Output: 10-14 mV Response 90%: 13 sec Expected Life: 60 months Warranty: 18 months	Range: 0-100% Oxygen Signal Output: Response 90%: 13 sec Expected Life: 60 months Warranty: 18 months
<i>ANI-75-KE7</i>	<i>ANI-75-KE8</i>	<i>ANI-75-KE9</i>	<i>ANI-77</i>
			
Range: 0-100% Oxygen Signal Output: Response 90%: 13 sec Expected Life: 60 months Warranty: 18 months	Range: 0-100% Oxygen Signal Output: Response 90%: 6 sec Expected Life: 60 months Warranty: 18 months	Range: 0-100% Oxygen Signal Output: Response 90%: 6 sec Expected Life: 60 months Warranty: 18 months	Range: 0-100% Oxygen Signal Output: 23-30 mV Response 90%: 6 sec Expected Life: 60 months Warranty: 18 months
<i>ANI-915</i>	<i>ANI-915-1</i>	<i>ANI-915-2</i>	<i>ANI-915-2i</i>
			
Range: 0-100% Oxygen Signal Output: 53-67 μ A Response 90%: 20 sec Expected Life: 60 months Warranty: 18 months	Range: 0-100% Oxygen Signal Output: 53-67 μ A Response 90%: 20 sec Expected Life: 60 months Warranty: 18 month	Range: 0-100% Oxygen Signal Output: 23-30 μ A Response 90%: 20 sec Expected Life: 60 months Warranty: 18 months	Range: 0-100% Oxygen Signal Output: Response 90%: 20 sec Expected Life: Warranty: 18 months

<p align="center"><i>ANI-915-3</i></p>	<p align="center"><i>ANI-915-4</i></p>	<p align="center"><i>ANI-915-G</i></p>	<p align="center"><i>ANI-917-J</i></p>
			
<p>Range: 0-100% Oxygen Signal Output: 53-67 μA Response 90%: 20 sec Expected Life: 60 months Warranty: 18 months</p>	<p>Range: 0-100% Oxygen Signal Output: 53-67 μA Response 90%: 20 sec Expected Life: 60 months Warranty: 18 months</p>	<p>Range: 0-100% Oxygen Signal Output: Response 90%: 20 sec Expected Life: 60 months Warranty: 14 months</p>	<p>Range: 0-100% Oxygen Signal Output: 7-13 mV Response 90%: 20 sec Expected Life: 60 months Warranty: 18 months</p>
<p align="center"><i>ANI-917-J1</i></p>	<p align="center"><i>ANI-917-J2</i></p>	<p align="center"><i>ANI-917-J3</i></p>	<p align="center"><i>ANI-917-M</i></p>
			
<p>Range: 0-100% Oxygen Signal Output: 7-13 mV Response 90%: 13 sec Expected Life: 60 months Warranty: 18 months</p>	<p>Range: 0-100% Oxygen Signal Output: 7-13 μA Response 90%: 13 sec Expected Life: 60 months Warranty: 18 months</p>	<p>Range: 0-100% Oxygen Signal Output: 7-13 mV Response 90%: 20 sec Expected Life: 60 months Warranty: 18 months</p>	<p>Range: 0-100% Oxygen Signal Output: 7-13 mV Response 90%: 13 sec Expected Life: 60 months Warranty: 18 months</p>
<p align="center"><i>ANI-917-M1</i></p>	<p align="center"><i>ANI-917-M2</i></p>	<p align="center"><i>ANI-917-MH</i></p>	
			
<p>Range: 0-100% Oxygen Signal Output: 7-13 mV Response 90%: 13 sec Expected Life: 60 months Warranty: 18 months</p>	<p>Range: 0-100% Oxygen Signal Output: Response 90%: 13 sec Expected Life: Warranty: 18 months</p>	<p>Range: 0-100% Oxygen Signal Output: 7-13 mV Response 90%: 13 sec Expected Life: 60 months Warranty: 18 months</p>	

<i>Brand</i>	<i>OEM Equipment</i>	<i>OEM Sensor Part #</i>	<i>Flotec Sensor Part #</i>
Air Shields (Dräger)	C200 / C2000	8362030	ANI-917-M
Analytical Industries	2000 Analyzer	AII-11-60	ANI-60
Analytical Industries	2000A Analyzer	AII-11-60	ANI-60
Analytical Industries	2000HC Analyzer	AII-11-60HC	ANI-60HC
Analytical Industries	2000M Analyzer	AII-11-60	ANI-60
Atom Medical	V-2200 Incubator		ANI917-MH1
Bird Products	6400	66030	ANI-917-J
Bird Products	8400	66030	ANI-917-J
Bio-med Devices		65004	ANI-33
Carefusion	SiPap Infant Flow Ventilator		ANI-917-MH
Carefusion	Vela		ANI-75-KE4
Criticare	Poet Monitor		ANI-917-MH
Criticare Systems	Poet	CAT644	ANI-917-M
Datascope	Anestar (new style)		ANI-917-J3
Datascope	Anestar (old style)	0600-00-0070	ANI-917-J2
Datascope	Monitor	0600-00-0002	ANI-917-M1
Dräger	Baby Log Vent		ANI-915-2i
Dräger	Cato	6850645	ANI-915-2
Dräger	C2000 Incubator		ANI-917-MH1
Dräger	Evita 2 & 4 Vent	6850645	ANI-915-2
Dräger	Evita XL		ANI-915-2i
Dräger	Fabius w/cable		ANI-915-1
Dräger	Fabius w/cable	6850645	ANI-915-2
Dräger	GS	6850645	ANI-915-2
Dräger	I8000 Incubator		ANI-915-2i
Dräger	Narkomed 1A		ANI-915-1
Dräger	Narkomed 1A	6850645	ANI-915-2
Dräger	Narkomed 2		ANI-915-1
Dräger	Narkomed 2	6850645	ANI-915-2
Dräger	Narkomed 2B		ANI-915-1
Dräger	Narkomed 2B	6850645	ANI-915-2
Dräger	Narkomed 3		ANI-915-1
Dräger	Narkomed 3	6850645	ANI-915-2
Dräger	Narkomed 4		ANI-915-1
Dräger	Narkomed 4	6850645	ANI-915-2
Dräger	Narkomed 6000		ANI-915-1
Dräger	Narkomed 6000	6850645	ANI-915-2
EME	Infant Flow CPAP		ANI-917-M
Evita	2 Vent		ANI-915-2i
Evita	4 Vent		ANI-915-2i
Evita	XL		ANI-915-2i

<i>Brand</i>	<i>OEM Equipment</i>	<i>OEM Sensor Part #</i>	<i>Flotec Sensor Part #</i>
Hamilton	Aladdin	51250	ANI-33
Hamilton	Amadeus	51250	ANI-33
Hamilton	Arabella CPAP	51254	ANI-917-MH
Hamilton	Galileo Ventilator	51254	ANI-917-MH
Hamilton	Ralphael Ventilator	51254	ANI-917-MH
Hamilton	Veolar	51250	ANI-33
Hudson	5525	5565	ANI-33-1
Hudson	5525	5565	ANI-33-2
Hudson	5550	5565	ANI-33-1
Hudson	5550	5565	ANI-33-2
Hudson	5552	5565	ANI-33-1
Hudson	5552	5565	ANI-33-2
Hudson	5577	5500	ANI-55
Hudson	5590	5500	ANI-55
Hudson	5828	5565	ANI-33-1
Hudson	5828	5565	ANI-33-2
Hudson	6477	5500	ANI-55
Hudson	5800	5803	ANI-58
Hudson	5801	5803	ANI-58
Hudson	5802		ANI-58HC
Hudson	RCI 5803	5804	ANI-58HC
Infrasonics	Adult Star Ventilator	60116	ANI-917-M
Lifecare (Respironics)	PLV-102	37-079	ANI-917-J
Maxtec	MaxBlend		ANI-75-KE9
Maxtec	OM-25A Analyzer		ANI-75-KE1
Maxtec	OM-25A Analyzer		ANI-75-KE9
Maxtec	OM-25A Analyzer		ANI-917-J1
Maxtec	OM-25E Monitor		ANI-917-J1
Maxtec	OM-25ME Monitor		ANI-75-KE9
Maxtec	OM-25MEL Monitor		ANI-75-KE9
Maxtec	OM-25MEL Monitor		ANI-917-J1
Maxtec	OMO2+A Analyzer		ANI-75-KE10+
Maxtec (Ceramatec)	OM25AE	MAX25E/250E	ANI-917-J1
MSA	Mini Ox I	406931	ANI-917-J
MSA	Mini Ox II	406931	ANI-917-J
MSA	Mini Ox 3	406931	ANI-917-J
MSA	Mini Ox 3000	406931	ANI-917-J
Newport	E360 Vent		ANI-75-KE7
Newport	E500 Vent		ANI-75-KE7
Newport	Wave Ventilator	OM100	ANI-77
Newport	Wave Ventilator	OM105	ANI-77
Ohmeda	4700 Oxicap	0237-2034-700	ANI-915

<i>Brand</i>	<i>OEM Equipment</i>	<i>OEM Sensor Part #</i>	<i>Flotec Sensor Part #</i>
Ohmeda	5100 Monitor	0237-2034-700	ANI-915
Ohmeda	5100 Monitor		ANI-915-3
Ohmeda	5120 Monitor	0237-2034-700	ANI-915
Ohmeda	5120 Monitor		ANI-915-3
Ohmeda	5140 Monitor	0237-2034-700	ANI-915
Ohmeda	5140 Monitor		ANI-915-3
Ohmeda	5150 Monitor		ANI-915
Ohmeda	5150 Monitor		ANI-915-3
Ohmeda	5250	0237-2034-700	ANI-915
Ohmeda	5250 RGM	0237-2034-700	ANI-915
Ohmeda	7100 Smart Vent	6050-004-110	ANI-915-4
Ohmeda	7200 Monitor	0237-2034-700	ANI-915
Ohmeda	7200 Monitor		ANI-915-3
Ohmeda	7800 Monitor	0237-2034-700	ANI-915
Ohmeda	7800 Monitor		ANI-915-3
Ohmeda	7810 Monitor	0237-2034-700	ANI-915
Ohmeda	7810 Monitor		ANI-915
Ohmeda	7900 Smart Vent	6050-004-110	ANI-915-4
Ohmeda	ADU w/OEM Holder	6050-004-110	ANI-915-4
Ohmeda	Aestiva 3000 / 5	6050-004-110	ANI-915-4
Ohmeda	Asys w/OEM Holder	6050-004-110	ANI-915-4
Ohmeda	Avance	6050-004-110	ANI-915-4
Ohmeda	Excel Series	0237-2034-700	ANI-915
Ohmeda	Giraffe Incubator		ANI-915-G
Ohmeda	Modules Series	0237-2034-700	ANI-915
Penlon	AV900		ANI-60-3
Penlon	Model AVS		ANI-60-3
Puritan Bennett	740 Ventilator (older version)	G-062010	ANI-75-KE2
Puritan Bennett	760 Ventilator	4-072214-00	ANI-917-M2
Puritan Bennett	840 Ventilator	4-072214-00	ANI-917-M2
Puritan Bennett	7820 Monitor	4-020933	ANI-917-J
Respironics	Esprit	1001454	ANI-917-J
Sechrist	650	43170	ANI-77
Sensidyne	1000E Analyzer		ANI-60-08
Sensidyne	5000M Monitor		ANI-60-23
Siemens	300 series Ventilator	6419332	ANI-917-M
Siemens	900 series Ventilator	9004979	ANI-917-M
SLE	2000 Ventilator		ANI-75-KE6
SLE	5000 Ventilator		ANI-75-KE8
SLE	T-801 monitor		ANI-917-J
Spacelabs	Gas Analyzer	90514	ANI-917-M
Stephan Christina	Ventilator	110760009	ANI-917-J

<i>Brand</i>	<i>OEM Equipment</i>	<i>OEM Sensor Part #</i>	<i>Flotec Sensor Part #</i>
Taema	KV050501		ANI-33
Taema (France)	Alys		ANI-917-J3
Taema (France)	Caesar (newer version)		ANI-917-J
Taema (France)	Caesar (older version)		ANI-917-M
Taema (France)	Clarys		ANI-917-J3
Taema (France)	Eolo	R17 (C-43690)	ANI-917-J
Teledyne	60T Analyzer	T-7 (A51327)	ANI-77
Teledyne	191 Monitor	T-7 (A51327)	ANI-77
Teledyne	200 Monitor	T-7 (A51327)	ANI-77
Teledyne	AX 300 Analyzer		ANI-917-J
Teledyne	MX 300 Monitor		ANI-917-J
Versamed	iVent		ANI917-M
Viasys	Avea Vent	62689	ANI-917-M
Viasys	Sentry Blender		ANI-75-KE9
Viasys	Vela	16101	ANI-75-KE4

<i>Replacement Brand</i>	<i>Replacement Sensor Part #</i>	<i>Flotec Sensor Part #</i>
Air Shields (Dräger)	8362030	ANI-917-M
Bio-Med Devices	65004	ANI-33
Bio-Tek	74233	ANI-33
City Technology	MOX-2	ANI-917-J
City Technology	MOX-4	ANI-77
City Technology	MOX-6	ANI-915-2
City Technology	MOX-9	ANI-917-M
Criticare Systems (Poet)	CAT644	ANI-917-M
EME	M672P	ANI-917-M
EnviteC	OOM101	ANI-915
EnviteC	OOM102/ OOM102-Air	ANI-917-M
EnviteC	OOM102-1	ANI-917-J
EnviteC	OOM103	ANI-917-M
EnviteC	OOM104	ANI-915-2
EnviteC	OOM105	ANI-77
EnviteC	OOM106	ANI-917-M1
EnviteC	OOM107	ANI-33
EnviteC	OOM107-2	ANI-33-2
EnviteC	OOM110	ANI-915-4
EnviteC	OOM201	ANI-915-2
EnviteC	OOM202	ANI-917-M
Fraser Lake	66-0453	ANI-33
Hamilton	51250	ANI-33
Hewlett Packard	15201-A	ANI-33
Hudson RCI	5557	ANI-915-2
Hudson RCI	5558	ANI-77
Hudson RCI	5565	ANI-33-2

<i>Replacement Brand</i>	<i>Replacement Sensor Part #</i>	<i>Flotec Sensor Part #</i>
Hudson RCI	5566	ANI-917-M
Hudson RCI	5567	ANI-917-J
Hudson RCI	5568	ANI-33
Infrasonics	60116	ANI-917-M
International Tech (Gambert)	M-01	ANI-915
International Tech (Gambert)	M-02 / MO8	ANI-915-2
International Tech (Gambert)	M-03	ANI-917-J
International Tech (Gambert)	M-04/M-12	ANI-917-M
International Tech (Gambert)	M-05	ANI-33
International Tech (Gambert)	M-06	ANI-33-2
International Tech (Gambert)	M-07	ANI-77
International Tech (Gambert)	M-10	ANI-915-4
International Tech (Gambert)	M-13	ANI-55
International Tech (Gambert)	M-25	ANI-917-J1
Libra	OM500	ANI-33
Lifecare (Respironics)	37-079	ANI-917-J
Maquet (Siemens)	6419332/9004979	ANI-917-M
Maxtec (Ceramatec)	MAX-1	ANI-915
Maxtec (Ceramatec)	MAX-2	ANI-915-1
Maxtec (Ceramatec)	MAX-8	ANI-915-3
Maxtec (Ceramatec)	MAX-9	ANI-55
Maxtec (Ceramatec)	MAX-10	ANI-915-4
Maxtec (Ceramatec)	MAX-11	ANI-915-2
Maxtec (Ceramatec)	MAX-12	ANI-917-M
Maxtec (Ceramatec)	MAX-13/ MAX-13-250	ANI-917-J
Maxtec (Ceramatec)	MAX-15	ANI-33-2
Maxtec (Ceramatec)	MAX-16	ANI-917-M
Maxtec (Ceramatec)	MAX-17/ MAX-17-250	ANI-77
Maxtec (Ceramatec)	MAX-18	ANI-33-1
Maxtec (Ceramatec)	MAX-19	ANI-33
Maxtec (Ceramatec)	Max-22C	ANI-917-M
Maxtec (Ceramatec) OM2	MAX-25/250	ANI-75-KE1
MSA (Catalyst Research)	406931	ANI-917-J
MSA (Catalyst Research)	478841	ANI-33
MSA (Catalyst Research)	485905	ANI-33-1
MSA (Catalyst Research)	655262	ANI-33-2
MSA (Catalyst Research)	655263	ANI-915
MSA (Catalyst Research)	655264	ANI-915-2
MSA (Catalyst Research)	655265	ANI-55
MSA (Catalyst Research)	716044 (PB 840)	ANI-917-M
MSA (Catalyst Research)	804674	ANI-77
Newport Medical	OM100/OM105	ANI-77
Pacifitech	PT-01	ANI-33
Pacifitech	PT-02	ANI-33-2
Pacifitech	PT-03	ANI-77
Pacifitech	PT-05	ANI-55
Pacifitech	PT-11/PT-11A	ANI-917-J

<i>Replacement Brand</i>	<i>Replacement Sensor Part #</i>	<i>Flotec Sensor Part #</i>
Pacifitech	PT-12A/ PT-12AIR	ANI-917-M
Pacifitech	PT-14	ANI-915
Pacifitech	PT-14A	ANI-915-2
Pacifitech	PT-20A	ANI-915-4
Pacifitech	PT-24A	ANI-915-1
Pacifitech	PT-25A	ANI-917-J1
Pacifitech	PT-28A	ANI-915-3
Pacifitech	PT-250	ANI-75-KE1
Pacifitech	PT-65	ANI-33-2
Paragon Services	PS-1000	ANI-915
Paragon Services	PS-2000	ANI-915-2
Physio Control	803000-18	ANI-915-1
Resprionics (Esprit)	1001454	ANI-917-J
Sechrist	43170	ANI-77
Sensor Tec (BPR Medical)	ST-01	ANI-33
Sensor Tec (BPR Medical)	ST-02	ANI-33-1
Sensor Tec (BPR Medical)	ST-03	ANI-77
Sensor Tec (BPR Medical)	ST-05	ANI-55
Sensor Tec (BPR Medical)	ST-10	ANI-915
Sensor Tec (BPR Medical)	ST-11	ANI-917-J
Sensor Tec (BPR Medical)	ST-12	ANI-917-M
Sensor Tec (BPR Medical)	ST-14	ANI-915-2
SLE T-801 monitor		ANI-917-J
Spacelab's	90514	ANI-917-M
Stephan Christina Ventilator	110760009	ANI-917-J
Taema (France)	YR010100	ANI-917-J
Taema (France)	YRO49700	ANI-917-M
Teledyne	C1R (B-39022)	ANI-33-1
Teledyne	C2R (B-39023)	ANI-33-2
Teledyne	R13 (B39636)	ANI-55
Teledyne	R17 (C-43690)	ANI-917-J
Teledyne	R22 (C-44611)	ANI-917-M
Teledyne	Class R-30	ANI-915-4
Teledyne	Class R-31	ANI-915-1
Teledyne	Class R-32	ANI-915-3
Teledyne	R15 (C-41131)	ANI-915
Teledyne	R23 (C-61339)	ANI-915-2
Teledyne	R24 (C-61340)	ANI-917-M
Teledyne	T-7 (A51327)	ANI-77
Teledyne	T1 (A-37016)	ANI-33
Vascular Technology	100580	ANI-917-J
Vascular Technology	101272	ANI-915
Vascular Technology	101372	ANI-915-2
Vascular Technology	103100	ANI-917-M
Viasys Avea	62689	ANI-917-M