

## Regulator, InGage



The Flotec **InGage™ Regulator** is a medical gas regulator with a unique internal pressure gauge built into the body of the regulator with a window to show the cylinder pressure. The gauge features a color coded 0-3000 PSI range that can be reset in the field. The all brass, high pressure, yoke connector has a 10-micron sintered brass inlet filter and a permanent o-ring type yoke washer seal. The regulator is a single stage, piston type, backpressure compensated, pressure reducer with an integral relief valve and a “click style” flowmeter that has 12 flow positions that will operate in any orientation or during extreme vibration. The unique flow rotor assembly has a second integrated 20-micron, sintered stainless steel filter built into the laser drilled orifice plate providing an absolute level of contamination protection. Currently available in nine different Medical Gases, 16 colors, six pressure settings, five International standards, two different metals and a nearly unlimited number of outlet connections. The **InGage™ Regulator** design allows the user to customize the regulator so that it will fit almost any requirement. The unique and distinctive Flotec Dodecagon (12 sided) shape helps prevent confusion with other brands, prevents the regulator from rolling off tabletops, and makes it easier to crack pecans, kill bugs and drive nails.

- Unique Internal Pressure Gauge
- Unique Field Adjustable Pressure Gauge
- Unique Laser Drilled Orifice Plate with Integral Filter
- Unique Dome Handle
- Unique O-Ring type Yoke Washer Seal
- Unique Swivel Yoke & Body
- Unique Glow in the Dark Flow Numbers
- Unique Dodecagon (12 sided) shape
- Unique Double Filter System
- Unique Thread Insert Design
- Unique Custom Laser Etching available
- MRI Compatible and Safe
- 9 Medical Gas configurations available in 13 Different Anodized Colors
- All Metal Construction (available in Anodized Aluminum and Plated Brass)
- FDA CGMP (CFR 21, Parts 803-820), ISO9001:2008, ISO13485:2003, and CE 0086
- **Manufactured exclusively in the United States of America**



Part Number	Description	List Price
DR8033-7003P4	CGA 870, 0-25 LPM, Hose Barb	174.00
DR8133-7003P4	CGA 870, One DISS 1240 Check Valve, 0-25 LPM, Hose Barb	184.00
DR8333-7003P4	CGA 870, Two DISS 1240 Check Valves, 0-25 LPM, Hose Barb	194.00
DR5103-6003	CGA 540, One DISS Check Valve, Pressure Reducer, Plain Cap	111.00
DR5303-6003	CGA 540, Two DISS Check Valves, Pressure Reducer, Plain Cap	121.00
DR5303-7003P4	CGA 540, Two DISS 1240 Check Valves, 0-25 LPM, Hose Barb	187.00

Match Ambulance outlets with Regulator outlets for quicker patient transfers

- o Ohio Quick Disconnect and DISS-1240 Check Valve
- o Chemetron Quick Disconnect and DISS-1240 Check Valve
- o DISS-1240 Check Valve



Part Number	Description	List Price
DR8E33-7003P4	CGA 870, One Ohio & One Check Valve, 0-25 LPM, Hose Barb	271.00
DR8Q33-7003P4	CGA 870, One Chemetron & One Check Valve, 0-25 LPM, Hose Barb	262.00

### Anodized Aluminum

Part Number	Description	List Price
DR8032-7002P4	CGA 870, 0-25 LPM, Hose Barb, Silver (Clear)	178.00
DR8033-7003P4	CGA 870, 0-25 LPM, Hose Barb, Green	174.00
DR8034-7004P4	CGA 870, 0-25 LPM, Hose Barb, Yellow	192.00
DR8035-7005P4	CGA 870, 0-25 LPM, Hose Barb, Blue	178.00
DR8036-7006P4	CGA 870, 0-25 LPM, Hose Barb, Red	182.00
DR8037-7007P4	CGA 870, 0-25 LPM, Hose Barb, Brown	192.00
DR8038-7008P4	CGA 870, 0-25 LPM, Hose Barb, Gold	196.00
DR8039-7009P4	CGA 870, 0-25 LPM, Hose Barb, Black	192.00
DR803A-700AP4	CGA 870, 0-25 LPM, Hose Barb, Purple	182.00
DR803B-700BP4	CGA 870, 0-25 LPM, Hose Barb, Orange	196.00
DR803C-700CP4	CGA 870, 0-25 LPM, Hose Barb, Camouflage Desert	212.00
DR803P-700PP4	CGA 870, 0-25 LPM, Hose Barb, Pink	182.00
DR803W-700WP4	CGA 870, 0-25 LPM, Hose Barb, Camouflage Woodland	212.00

### Plated Brass

Part Number	Description	List Price
DR803T-700TP4	CGA 870, 0-25 LPM, Hose Barb, Silver (Tri-Metal)	230.00
DR8033-7003P4	CGA 870, 0-25 LPM, Hose Barb, 23k Gold Plated	240.00

