

## Minireg Series Oxygen Cylinder Regulators

Five Year Warranty

MRI Compatible and Safe

ML800-280xx



MK800-280xx



MK800- 663x



ML500-280xx



MK500-280xx



MK500-663x



All Flotec Medical Gas Regulators feature 10 micron high-pressure bronze sintered inlet filtration and an all brass high pressure internal gas pathway.

The flow meter module offers 12 flow positions and unlike Thorpe tube styles operates accurately in any position or while moving.

Flow patency and accuracy are maintained with our unique 20 micron "Rotor Filter" positioned proximal to the flow orifices.

Available in CGA 540 and CGA 870

Custom Laser Engraving available on all regulators

Heavy Duty Service with outstanding impact resistance

All-Metal-Construction with unibody window design

True Finger-Tight Post Valve Seal

Anti-Blowout Design utilizing no snap rings

Flow ranges: 0-6 LPM, 0-8 LPM, 0-15 LPM, 0-25 LPM

Low Profile Gauge with impact resistant rubber boot

Large view window with high visibility Night Glow Decals

**Dome Style Knob  
Handle Available**



**on ALL Flotec  
Yoke Body Regulators!**