

M L 8 0 0 5 -- 2 8 0 5 M2

B-Regulator, Mini Code (2013-04-14).123.pdf

FLOW RANGE:
See B-Flows
chart for part
number and
list price

COLOR OPTION:
NONE - Green Anodized Body
5 - Blue Anodized Body

-OR-
FITTING OUTLET:
A - 1 LPM Orifice
B - 2 LPM Orifice
C - 3 LPM Orifice
D - 4 LPM Orifice
E - 5 LPM Orifice
F - 6 LPM Orifice
H - 8 LPM Orifice
J - 10 LPM Orifice
L - 12 LPM Orifice
M - 15 LPM Orifice
N - 20 LPM Orifice
P - 25 LPM Orifice
T - 60 LPM Orifice
V - 70 LPM Orifice

FLOW CARTRIDGE OPTIONS:
0 - No Flow Cartridge Options
A - 20 PSI Outlet Pressure

-OR-
CAP OPTIONS:
0 - Plain Cap, No Outlets
1 - Hose Barb 1/8" Pipe
2 - 1/8" Female Pipe Port
P - Short Cap w/ 1/8" FNPT

WINDOW & CONTROLLED FLOW OUTLET POSITIONS:

- 0 - Window @ 0 & Hose Barb @ 180
- 1 - Window @ 0 & DISS-1240 fitting @ 180
- 8 - Window @ 90 & Hose Barb @ 180
- 9 - Window @ 90 & DISS-1240 fitting @ 180

FLOW CARTRIDGE CONFIGURATION:
2 - Large Knurled Knob
6 - Cap

-OR-
CAP OUTLET POSITIONS:
0 - Cap Outlet Position on Centerline
2 - Cap Outlet Position @ 0 Radial
3 - Cap Outlet Position @ 90 Radial
6 - Cap Outlet Position @ 180 Radial
9 - Cap Outlet Position @ 270 Radial
D - Cap Outlet Position @ 0, 90, 180, 270 Radial

COLOR OPTION:
NONE - Green Anodized Body
5 - Blue Anodized Body

BODY OPTION(S):
0 - No Body Options
3 - Dome-Handle
5 - Dome-Handle & 20 PSI Outlet Pressure
A - 20 PSI Outlet Pressure

PRESSURE REDUCED BODY OUTLET(S):
0 - No Pressure Body Outlets
1 - One DISS-1240 Check Valve
2 - One 1/8" FNPT (Plugged)
3 - Two DISS-1240 Check Valves

INLET CONNECTION:
5 - CGA 540 Nut & Nipple
8 - CGA 870 Yoke

GAUGE & BODY CONFIGURATIONS:
K - Left Hand, No Gauge
L - Left Hand, 1.4", 0 to 3000 PSI Gauge @ 60 degrees
R - Right Hand, 1.4", 0 to 3000 PSI Gauge @ 300 degrees

PRODUCT GROUP:
M - Mini Regulator