

Reconditioned InGage Regulators

Flotec, Inc. has limited supply of InGage regulators available at a special price. These regulators were part of a completed, long-term rental program.



This regulator, DC8B03-7003P1 when purchased new, listed for \$225.00, and is currently available fully reconditioned, for the low price of **\$ 59.00** with a two-year warranty. All the filters, seals and orifices have been replaced. This regulator features a Small T-Handle, a single DISS-1240 Check Valve, and Hose Barb, a 0 to 25 flow range (0, 0.25, 0.5, 1, 2, 3, 4, 6, 8, 10, 15, 25), with a Large Diameter Knob.

- ✓ **Internal Pressure Contents Gauge:** No external gauge to break on impact
- ✓ **Full 0-3000 psi Pressure Range:** Composite and HP Aluminum Cylinder Ready
- ✓ **Easy to read Tri-Color Coded Pressure Range Indicators.**
- ✓ **Precision flow rates using optically laser drilled orifices.**
- ✓ **All-Metal-Construction assures long service life.**
- ✓ **Permanent Laser engraved markings**

