

Reconditioned InGage Regulators

Flotec, Inc. has limited supply of InGage regulators available at a very special price. These regulators were part of a long term rental program that has ended.



This regulator, purchased new, lists for \$225.00, but is currently available, reconditioned, for the low price of \$79.00 with a T-handle or \$89.00 with a dome handle (as pictured). Available in **Pink** and **Green**.

This new exclusive “**Pink**” regulator, only available through Flotec, comes with a CGA 870 connection, DCV-1240, hose barb outlet, 0-25 LPM, with reconditioned components. It carries a 5 year warranty.

- ✓ **Internal Pressure Contents Gauge:** No external gauge to break on impact
- ✓ **Full 0-3000 psi Pressure Range:** Composite and HP Aluminum Cylinder Ready
- ✓ **Easy to read Tri-Color Coded Pressure Range Indicators:** Visible from as far as 8 to 10 feet away.
- ✓ **Precision flow rates using optically laser drilled orifices.**
- ✓ **Finger-Tight Dome Handle CGA870 Inlet Connection. No wrenching!**
- ✓ **All-Metal-Construction assures long service life.**
- ✓ **Permanent Laser engraved markings**
- ✓ **MRI Compatible and Safe.**



Flotec, Inc. Order Form

Date: _____

Person Placing Order: _____

Company: _____

Ship- to Address: _____

Phone: _____ Fax: _____

e-mail: _____

Would you like your confirmation sent by fax or email? _____

| Part # | Qty. | Price | Extended | |
|---------------------------------------|-------|----------|----------|--------------|
| DC8B03-7003P5 Green with T-Handle | _____ | \$ 79.00 | _____ | |
| DC8B33-700P3P5 Green with Dome-Handle | _____ | \$ 89.00 | _____ | |
| DC8B0P-700PP5 Pink with T-Handle | _____ | \$ 79.00 | _____ | |
| DC8B3P-700PP5 Pink with Dome-Handle | _____ | \$ 89.00 | _____ | |
| Freight: Pre paid and added | | | _____ | TOTAL |

Name on Card: _____

Address (Bill To) _____

| Card Type | Exp Date | V Code | Card Number |
|-----------|----------|--------|-------------|
| | | | |

Please fax your order form to (800) 515-9254 or call your order in to our order desk at (800) 401-1723.
A confirmation will be sent back to you via your preference of fax or email.