





## Ambulance Respiratory Systems

Flotec has developed a new Respiratory System designed specifically for Ambulance Applications. The features and functions of the Regulator and the Flow Control Valve are matched for increased efficiency.

- Regulator and Flow Control Valve have matching flow ranges covering the full range of 0-25 LPM.
- Regulator and Flow Control Valve have matching hose barb outlets.
- Regulator and Flow Control Valve have optional matching DISS-1240 Check Valve auxiliary outlets.
- Flowmeter has an optional removeable Resident Nipple Adapter for CPAP or Humidifier applications.
- Flow Control Valve has multiple outlet options to match Regulator.



DR8133-7003P4




FN0-2EEP6

FN0-2EOP6



T10-BEED3P6

- Dual Flow Control Valve from a single wall outlet.
- Optional Pediatric plus High Flow Control Valve.
- Optional removeable Resident Nipple Adapter for CPAP or Humidifier applications.
- Optional DISS-1240 Check Valve, Female Ohio or Chemetron.
- Optional T or Y styles.





### Ambulance Set-Up 0-25 LPM (no DISS check valve)

| Part Number   | Description   |
|---------------|---|
| F10-2E0P6     | Flow Control Valve 0-25 LPM, DISS Outlet              |
| FN0-2E0P6     | Flow Control Valve 0-25 LPM, DISS Outlet with HBD2    |
| T10-BE0D3P6   | Dual Flow Control Valve 0-4 & 0-25 LPM                |
| DR8033-7003P4 | CGA 840, Regulator 0-25 LPM                           |
| DR5103-6003   | CGA 540, Pressure Reducer, 1 DISS CV                  |
| RZ510-600     | CGA 540, Pressure Reducer, Transducer port, 1 DISS CV |



DR8033-7003P4



FN0-2E0P6

### Ambulance Set-Up 0-25 LPM (1 DISS check valve)

| Part Number   | Description   |
|---------------|---|
| F10-2EEP6     | Flow Control Valve 0-25 LPM, DISS Outlet                      |
| FN0-2EEP6     | Flow Control Valve 0-25 LPM, 1 DISS CV, DISS Outlet with HBD2 |
| T10-BEED3P6   | Dual Flow Control Valve 0-4 & 0-25 LPM, 1 DISS CV             |
| DR8133-7003P4 | CGA 840, Regulator 0-25 LPM, 1 DISS CV                        |
| DR5103-6003   | CGA 540, Pressure Reducer, 1 DISS CV                          |
| RZ510-600     | CGA 540, Pressure Reducer, Transducer port, 1 DISS CV         |



DR8133-7003P4



FN0-2EEP6

### Ambulance Set-Up 0-25 LPM (2 DISS check valve)

| Part Number   | Description   |
|---------------|---|
| FN0-2EEP6     | Flow Control Valve 0-25 LPM, 1 DISS CV                        |
| FN0-2EEP6     | Flow Control Valve 0-25 LPM, 1 DISS CV, DISS Outlet with HBD2 |
| T10-BEED3P6   | Dual Flow Control Valve 0-4 & 0-25 LPM, 1 DISS CV             |
| DR8333-7003P4 | CGA 840, Regulator 0-25 LPM, 2 DISS CV                        |
| DR5303-6003   | CGA 540, Pressure Reducer, 2 DISS CV                          |
| RZ530-600     | CGA 540, Pressure Reducer, Transducer port, 2 DISS CV         |



DR8333-7003P4



FN0-2EEP6

| P/N   |  |
|-------|--|
| HBD24 |  |
| HBD26 |  |
| HBD2  |  |





## CGA 540 Pressure Reducers for M-Cylinders

- Hand Tight and Wrench styles.
- All styles include high pressure O-Ring.
- Light weight.
- One or Two DISS-1240 Check Valves.



DR5103-6003



DR5306-6009

## Transducer Port Configurations



RZ510-600



RZU10-600



RZU00-600E